

Mounting Consumables & Accessories



Pigments

Use with castable resins for color coding or creating contrast

- 20-8505 Black, 1.5oz [45mL]
- 20-8506 Red, 1.5oz [45mL]
- 20-8507 Blue, 1.5oz [45mL]



Release Agent

Liquid release agent for easier removal of mounts from castable molds or compression mounting presses

- 20-8186-004[†] 4oz [120mL]
- 20-8186-032[†] 32oz [950mL]



Mold Release Spray

Spray release agent for use on castable mounting molds

- 20-3050-008 8oz [0.24L]



Mold Release Powder

Powder release agent for use on mounting presses

- 20-3048 2oz [45g]



SamplKlip

Stainless Steel support clip for use with all mounting systems.*

- 0.25 H x 0.550 W x 0.350in L [6 x 14 x 9mm] 0.575g
- 20-4000-100 (Qty 100)



SamplKlip I

Plastic support clip best for castable mounting systems.*

- 0.25 H x 0.475 W x 0.3in L [~6 x 12 x 8mm] 0.230g
- 20-4100-100 (Qty 100)
- 0.25 H x 0.425 W x 0.25in L [~6 x 11 x 6mm] 0.230g
- 20-4100-100S (Qty 100)



Specimen Support Clip

Plastic support clip best for castable mount systems.†

- 0.25 H x 0.290 W x 0.375in L [6 x 7 x 9.5mm] 0.145g
- 20-4001-000 (Qty 1000)



UniClip Support Clip

Plastic support clip for use with all mounting systems.

- 0.4 H x 0.360 W x 0.500in L [10 x 9 x 13mm] 0.290g
- 20-5100-100 Clear (Qty 100)
- 113043 Black (Qty 100)



Plastic Mixing Cup

Graduated Plastic Mixing cup for mixing castable mounting systems

- 20-8176-100 8.5oz [250mL] (Qty 100)



Paper Mixing Cup

Non Graduated Paper Mixing cup for mixing castable mounting systems

- 20-8177-100 5oz [148mL] (Qty 100)



Stirring Sticks

For stirring castable mounting systems

- 20-8175 (Qty 1000)



MetKleer™ Adhesive Bases

For use with ring forms and castable systems

- 20-8188 4 x 5in [102 x 127mm] (Qty 10)



Conductive Filler

Nickel-based filler makes castable mounting systems conductive

- 20-8500 2 lb [0.9kg]



Flat Edge Filler

Enhances edge retention in castable mounting systems.

- 20-8196 1 lb [0.45kg]



Thermoplastic Cement

For adhering samples to glass slides or other specialty holders

- 40-8100 (Qty 12)



Crystalbond Mounting Wax

For adhering samples to glass slides or other specialty holders

- 40-8150 Flows at 165°F [74°C]
- 20-8145 Flows at 120°F [49°C]

* Compatible with specimens up to 0.200in [5mm] thick

† Compatible with specimens between 0.0035 - 0.090in [0.9 - 2.3mm]